

“DAIHAN” Cork Stopper, Ecologically Friendly Materials NEW

콜크마개, 친환경재료, Air Permeable, Superior Grade Corks

- **Cork Origin : 포르투갈(Portugal)**
- **0.5~2mm 콜크 입자를 사용하여 제작**
- **100% 천연Cork 재료를 사용하여 친환경적이며 인체에 무해함**
- **공기투과/통기성 “Air Permeable”**
- **재사용, 재활용 가능 “Renewable and Recyclable”**

Application & Features / 사용분야 및 특징

- * 완벽한 셀링 “Tight Fitting”
- * 단열 효과 “Insulation Effect”
- * 강한 신축성 “Strong Elastic”
- * 미생물, 화학, 농산물, 의학 및 기타 분야용
- * General Lab & Medical 용, Sample Vessel 및 기타 산업용
- * 콜크 강도 : 170~240 (kg/m³)
- * 가벼운 무게 “Light Weight”
- * Cell / Tissue-Culture 용
- * Bottle, Flask, Test-Tubes, 등의 마개용.



100% Natural Cork Stopper

Specifications & Ordering / 규격 및 주문

* The Size : Similar to Silicone Stopper Size

Cat. No.	Stopper Number	Stopper Size φ 직경 상 / 下 × 높이(mm)	Pk.	Price ₩
DH.ST802	2 호	Normal Size φ13 / 9 × h15mm	10ea	1,250
DH.ST803	3 호	φ15 / 11 × h20mm	10ea	1,380
DH.ST804	4 호	φ17 / 13 × h20mm	10ea	1,500
DH.ST804.1	4-1 호	φ19 / 15 × h20mm	10ea	1,630
DH.ST805	5 호	φ21 / 16 × h27mm	10ea	2,250
DH.ST806	6 호	φ23 / 18 × h27mm	10ea	2,500
DH.ST807	7 호	φ25 / 20 × h27mm	10ea	2,750
DH.ST808	8 호	φ27 / 21 × h30mm	10ea	3,250
DH.ST809	9 호	φ30 / 24 × h30mm	10ea	3,750
DH.ST811	11 호	φ33 / 27 × h30mm	10ea	4,250
DH.ST812	12 호	φ37 / 30 × h39mm	10ea	6,000
DH.ST814	14 호	φ41 / 33 × h39mm	5ea	3,380
DH.ST815	15 호	φ45 / 37 × h39mm	5ea	4,000
DH.ST816	16 호	φ50 / 42 × h39mm	5ea	4,750
DH.ST818	18 호	Large Size φ55 / 47 × h39mm	5ea	5,380
DH.ST820	20 호	φ60 / 52 × h39mm	5ea	6,130
DH.ST821	21 호	φ65 / 57 × h39mm	5ea	6,500
DH.ST823	23 호	φ70 / 62 × h39mm	1ea	1,500
DH.ST824	24 호	φ75 / 67 × h39mm	1ea	1,780



Normal Size(φ13~50)



Large Size(φ55~75)

Silicone Sleeve Septa / Stopper with Snap-lip

실리콘 슬리브 스냅 셉타 / 스톱퍼, 주사기 인젝션용 마개

Sleeve septa/ stopper with hollow plugs, rounded corners, and needle guide.
stopper folds down over neck of bottle. Is easily perforated with hypodermic needles.

DH.ST1146 Sleeve Septa, Rubber, 14/23, 14/35 공용
DH.ST1146A , 19/26, 14/38 공용
DH.ST1145 , 24/29, 24/40 공용

Pk.	Price, ₩
100ea	43,750
100ea	50,000
100ea	52,500



Silicone Sleeve Septa/Stopper

“Simport” Santoprene-Rubber Stoppers, for φ12~13mm & φ16mm tubes 산토프렌 러버 스톱퍼, Flexible, 냉동·냉장 보관 가능, Ideal for Needle Injection

SP.T407.12B Cap, for φ12~13mm tube, Blue
SP.T407.12Y , for φ12~13mm tube, Yellow
SP.T407.16B , for φ16mm tube, Blue
SP.T407.16Y , for φ16mm tube, Yellow

Pk.	Price, ₩	Cs.	Price, ₩
1000ea	96,930	10pk	940,221
1000ea	96,930	10pk	940,221
1000ea	97,830	10pk	948,951
1000ea	97,830	10pk	948,951



Santoprene-Rubber Stopper