




DURAN® The Highest Quality NMR Tubes NEW

NMR 튜브, Economic / Professional / Scientific 3 종류, High Quality

■ Supplied by DURAN® in Germany

- * NMR tubes are available, according to requirement, in three accuracy classes.
- * The correct tube can be selected depending on magnetic field and spin.
- * These tubes are noteworthy for their close tolerances and accuracy.
- * Quick and accurate test results are achievable.
- * DURAN® Borosilicate Glass α3.3 <특급 유리> * 완전 세척되어 있음.
- * 모든 Tube는 Laser check system을 거쳐 규격이 개별 테스트되며, 자성물질을 함유하지 않음.
- * Diameter와 Wall thickness가 일정하고 Straightness가 우수하며, Wilmad 제품과 동등이상의 품질.
- * Tube가 Thin wall이므로 극히 가는 Tube지만 많은 양의 Sample을 분석 가능.
- * 3가지 종류 : 품질에 따라 Economic / Professional / Scientific * Plug 포함.

주문 (Ordering)

(1) Economic NMR Tube

DU.231700117 NMR Tube, Economic type, OD $\phi 4.95 \pm 0.05\text{mm}$,
ID $\phi 4.20 \pm 0.05\text{mm}$, h 178mm, Camber 0.07 mm, 300 MHz

Ea. Price, ₩	Pk.	Price, ₩
3,220	50ea	155,390

(2) Professional NMR Tube

DU.231700211 NMR Tube, Professional type, OD $\phi 4.97 \pm 0.025\text{mm}$,
ID $\phi 4.20 \pm 0.025\text{mm}$, h 178mm, Camber 0.03 mm, 400 MHz

Ea. Price, ₩	Pk.	Price, ₩
13,440	25ea	325,620

(3) Scientific NMR Tube

DU.231700314 NMR Tube, Scientific type, OD $\phi 4.97 \pm 0.013\text{mm}$,
ID $\phi 4.20 \pm 0.025\text{mm}$, h 178mm, Camber 0.013 mm, 500 MHz

Ea. Price, ₩	Pk.	Price, ₩
24,990	5ea	121,180

Spare Caps for NMR Tubes

DU.2971701 Cap for NMR Tube, Blue
DU.2991702 , Red
DU.2991703 , Yellow

Ea. Price, ₩	Cs.	Price, ₩
1,000	250ea	206,990
1,000	250ea	206,990
1,050	250ea	215,990