

Cat. No	Description
---------	-------------

Articles

“Witeg” Colored Lab Bottles, Graduated in a Color-code system
칼라코드 래보라토리 바틀, Yellow & Red, 눈금부 & GL45 Screw Cap/Pour 링, ISO

		Ea.	Price, ₩
WI.5526100GE	Bottle, Yellow, 100 ml, GL45, $\phi 56 \times h100$ mm		12,480
WI.5526250GE	, 250 ml, GL45, $\phi 70 \times h138$ mm		13,500
WI.5526500GE	, 500 ml, GL45, $\phi 86 \times h176$ mm		17,250
WI.5526001GE	, 1000 ml, GL45, $\phi 101 \times h225$ mm		25,070
WI.5526002GE	, 2000 ml, GL45, $\phi 136 \times h260$ mm		55,440
WI.5526100RO	Bottle, Red, 100 ml, GL45, $\phi 56 \times h100$ mm		12,480
WI.5526250RO	, 250 ml, GL45, $\phi 70 \times h138$ mm		13,500
WI.5526500RO	, 500 ml, GL45, $\phi 86 \times h176$ mm		17,250
WI.5526001RO	, 1000 ml, GL45, $\phi 101 \times h225$ mm		25,070
WI.5526002RO	, 2000 ml, GL45, $\phi 136 \times h260$ mm		55,440



Colored Lab Bottles

“Wheaton” Media Lab Bottles, with Fine Graduation & White Marking Area
메디아 바틀, 눈금부 / 백색 Marking-area부, ASTM/EPA/USP

Graduated media bottles have a large opening, no-drip pour lip and sloping shoulders to facilitate pouring and cleaning. Bottles, rubber lined and Teflon lined caps are autoclavable.

- * 사용도 : Culture Media, Reagent, Sampling, Storage... 기타용도.
- * With Rubber lined Cap : Widely used for mixing and storage of culture media.
- * With Teflon lined Cap : Teflon liners are inert to most chemicals.

With Rubber Lined Phenolic Cap / 러버 라인드 캡 포함

		Ea. Price, ₩	Cs.	Price, ₩
WH.219755	Bottle, Media, 125 ml, $\phi 55 \times h123$ mm, Cap 33-430	5,680	48ea	252,000
WH.219757	, 250 ml, $\phi 67 \times h152$ mm, Cap 33-430	7,440	48ea	330,300
WH.219759	, 500 ml, $\phi 88 \times h192$ mm, Cap 33-430	11,140	24ea	247,480
WH.219760	, 1000 ml, $\phi 103 \times h235$ mm, Cap 38-430	14,980	24ea	332,530

With Teflon Lined Phenolic Cap / 테프론 라인드 캡 포함

		Ea. Price, ₩	Cs.	Price, ₩
WH.219815	Bottle, Media, 125ml, $\phi 55 \times h123$ mm, Cap 33-430	6,130	48ea	272,130
WH.219817	, 250ml, $\phi 67 \times h152$ mm, Cap 33-430	7,790	48ea	345,970
WH.219819	, 500ml, $\phi 88 \times h192$ mm, Cap 33-430	12,150	24ea	269,940
WH.219820	, 1000ml, $\phi 103 \times h235$ mm, Cap 38-430	15,720	24ea	348,990



Media Lab Bottles

“DAIHAN” DNA Hybridization Bottles, with Red PBT Cap
DNA 혼성화를 위한 튜브, with GL45 screwcap/PTFE liner, Boro 3.3 glass

Useful in nucleic acid and protein blotting procedures where a minimum of probe solution is utilized. These fit many hybridization ovens and a variety of bottle size are available to assure proper membrane fit.

- * Nucleic Acid, Protein Blotting에 사용.
- * 재질 : Borosilicate Glass $\alpha 3.3$
- * 적색 PBT cap with PTFE -protected seals
- * 대부분의 Hybridization Oven에서 사용가능.
- * Glass Wall Thickness : 2.3 mm

With GL45 Screwcap & Dia. $\phi 40$ mm, with Flat -bottom, Ideal for Hybridization Bottles, 90~ 280ml

		Ea. Price, ₩	Cs.	Price, ₩
DH.BO10011	Bottle, w/Red PBT GL45 cap, 90ml, $\phi 40 \times h100$	36,880	10ea	357,630
DH.BO10012	, 140ml, $\phi 40 \times h150$	36,880	10ea	357,630
DH.BO10013	, 190ml, $\phi 40 \times h200$	46,840	10ea	454,360
DH.BO10014	, 240ml, $\phi 40 \times h250$	47,350	10ea	459,290
DH.BO10015	, 280ml, $\phi 40 \times h300$	48,110	10ea	466,660




GL45 screwcap Bottles

“Wheaton” Glass Roller Culture Bottles, with Screwcaps
배양병, Roller Culture Apparatus용, Wide-neck id 40 mm, ASTM/USP

Bottles are made to exacting tolerances to provide higher cell yields. Glass are specially selected for optical clarity, interior surface smoothness and absence of glass imperfections. Wheaton Boro33 glass

- * Larger opening provides easier access to inside of vessel.
- * Screw Cap Size : 51-400(shallow skirt) with “SB” rubber liner

		Cs.	Price, ₩
WH.348522	Bottle, Roller, W-N, 1760 ml, $\phi 110 \times h270$ mm Surface Area : 840 cm ² , Working Area : 700 cm ²	4ea	455,140
WH.348524	Bottle, Roller, W-N, 2500 ml, $\phi 110 \times h355$ mm Surface Area : 1070 cm ² , Working Area : 940 cm ²	4ea	501,840
WH.348528	Bottle, Roller, W-N, 4230 ml, $\phi 110 \times h535$ mm Surface Area : 1680 cm ² , Working Area : 1555 cm ²	4ea	629,400



Wide Neck Culture Bottles