

Articles	Cat. No	Description
----------	---------	-------------

DAIHAN-brand® Multifunction Vortex Mixer Set, "VM-10", Up to 3,300rpm, with Certi. & Traceability with ϕ 76mm(3 inch) Platform Head & Pop-off Cup Head, Variable Speed Control, Continuous or Touch-on Mixing
다용도 볼텍스 믹서, 다양한 속도 조절, 강력하고 안정된 동작과 저소음



Multifunction Vortex Mixer Set "VM-10" with ϕ 76mm(3") Platform (PM110)

- **Made in Korea** ■ **CE Certified**
- **with Certificate & Traceability : Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System**
- **Ideal for Mixing Solutions in Test Tubes or Small Flasks with Circular, Oscillating Movement**
- **Variable Speed Control Allows Slow Speed Shaking Action up to High Speed Vortexing**


■ **PL (Product Liability) Insurance**

- * Gentle Mixing to re-suspending Pellet
- * Ideal for Mixing Solutions in Test Tubes or Small Flasks with Circular, Oscillating Movement
- * Supplied with ϕ 76mm Platform and Pop-Off Cup
- * Realize High Power, Stable Motion with Low Noise
- * Tough Metal Housing Provides a Stable Platform for all types of Mixing
- * 3-position Power Switch Operates with Continuous Action or "Touch On" Control
- * Variable Speed Control Allows Slow Speed Shaking Action up to High Speed Vortex
- * Aluminum Die-casting with Powder Coating
- * Various Shapes, Sizes and Materials of the Head allows for mixing of almost all Common Tubes or Containers

Specifications & Ordering information

Spec.	Models	"VM-10"
Speed Range		0 ~ 3,300 rpm
Shaking Motion		Orbital
Control		Analog Phase Controller
Head		included Pop-off Cup Head and ϕ 76mm Platform
Material		Aluminum Die-casting with Powder Coating
Dimensions (w × d × h) (mm)		160 × 128 × h145 mm
Weight		3 kg
Others		CE Certified 3-position Power Switch Operates Aluminum Die-casting with Powder Coating included Pop-off Cup Head and ϕ 76mm Platform
Packing Dimension (w × d × h) (mm)		260 × 200 × h220 mm
Gross Weight		4 kg
Power Consumption		30 W
Electrical Reqts.		1 Phase, AC 120V, 60Hz or AC 230V, 50 / 60Hz
Cat. No.	230V, 60Hz	DH.WVM00010
	230V, 50Hz	DH.WVM00015
	120V, 60Hz	DH.WVM00020
	120V, 50Hz	DH.WVM00025
Price, Each. ₩		158,000

* Other Specifications are available upon Customer's Request.



(1) Platform Head, ϕ 76mm(3inch) "PM110" (Included)

(2) Pop-off Cup Head "PM210" (Included)



(3) Microtube Platform Head, 0.5-2.0ml Tube "PM310" (Optional)

Remove Applicator (Included)

Spare Platform Heads

DH.WVM00110	(1) Platform Head, "PM110", ϕ 76mm(3inch)	Ea. Price, ₩ 11,000
DH.WVM00210	(2) Pop-off Cup Head, "PM210"	19,000

Optional Platform Heads

DH.WVM00310	(3) Microtube Platform Head(0.5~2.0ml), "PM310"	Ea. Price, ₩ 19,000
--------------------	---	-------------------------------