

Cat. No	Description	Articles
---------	-------------	----------

DAIHAN-brand® Digital Laboratory Freezer, “WLF”, -35℃ to +10℃, 420 Lit., with Certi. & Traceability with 3 × Shelves, Digital Fuzzy Control System, Fast Freezing Function, RS232C Interface, CFC-Free Refrigeration
디지털 저온 냉동고, 디지털 퍼지 컨트롤, 친환경 냉매(134A) 채용, 급속 냉동 기능

- with Certificate & Traceability : Controlled by Serial Number, Certificate, Delivery-information, and Traceable Data Base System
- Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- Fast Freezing Function
- Patented Jog-Shuttle Control System
- CFC-Free (134A) Refrigeration System
- CE Certified
- PL(Product Liability) Insurance

- * Suitable for Storage of Bio Reagent, Seed, Enzyme and Various Samples
- * Digital Fuzzy Control System Implementing Superior Temperature Accuracy
- * Storage Function : Set Values of Temperature and Timer
- * Digital LCD with Back-Light Function
- * Locking Mode supported for Experimental Safety (Input to Jog-Shuttle can be Disabled)
- * Intelligent Control of Refrigerator for Preventing Overload of Compressor
- * Corrosion Resistant 304 Stainless steel Chamber and Shelves
- * Over Current Protection, Sensor Error Detection, Leakage Breaking

Specification

Capacity & Model	420 Lit. Digital Laboratory Freezer “WLF-420”
Type	Upright Type
Temp. Range & Sensor	-35℃~+10℃, PT100
Temp. Resolution	1℃-Display, 1℃-Control
Refrigerant & Compressor	CFC-Free(134A) Refrigeration System, 1 HP
Display	Digital LCD with Back-Light Function
Controller	Digital Controller with Jog-Shuttle Switch (Turn+Push)
Included Shelf & Load per Shelf	3ea Stainless steel Shelves (# 304), 16 kg
Safety Device	Over Current Protector, Sensor Error Detector, Leakage Breaker
Material	Stainless steel (# 304)-Internal, Powder Coated steel-External
Dimension : Interior & Exterior	600×700×h1,000mm, 800×800×h1,700mm
Packing Size & Gross Weight	1,160×1,060×h2,160mm, 314 kg
Power Consumption	1.8 kW
Power Supply*	1 Phase, AC 120 V, 60 Hz or AC 230 V, 50 / 60 Hz

* Other Specifications are available upon Customer's Request.

※ 운임 약 ₩ 55,000 별도!

Ordering Information

DH.WLF00420	Lab Freezer, 230V, “WLF-420”	Ea. Price, ₩	6,302,900
DH.WLF01420	, 120V, “WLF-420”	Ea. Price, ₩	6,302,900
Spare Flat Shelf, Stainless steel			
DH.WLF10420	Flat Shelf, w580×d680mm, “LFS420”, for “WLF-420”	Ea. Price, ₩	82,640
Validation (IQ, OQ) Service			
DH.WLF30420	Validation Service(IQ, OQ), “LFV420”, for “WLF-420”	Ea. Price, ₩	문의요함



Digital Laboratory Freezer “WLF-420”



with Inner Doors “WLF-420”