

## Exceptional fit and feel for extended wear



superior tactile  
sensitivity

### FEATURES AND BENEFITS:

#### Ultimate in Tactile Sensitivity

Featuring excellent tactile sensitivity, trust the Microflex<sup>®</sup> 63-854 for precise handling of sensitive or delicate items, such as laboratory instruments.

#### Reliable and Confident Grip in Wet or Dry Conditions

Fully textured fingertips provide consistent and secure grip in messy work environments – perfect for either wet or dry applications.

#### Designed for Comfort and Productivity

Superior elasticity means gloves move with you for greater comfort while you work. Designed to be easy to slip on and remove – to help improve productivity.

#### Outstanding Quality

Exceeds standards for tensile strength, elasticity/elongation and other quality standards for average latex gloves. Consistent performance provides reliable protection across varied industries.

#### Industries

- Life Sciences (Pharmaceutical, Biopharmaceutical & Medical Devices)
- Electronics
- Laboratory

#### Applications

- Analytical testing/measurements
- Assembly and inspection
- Biotechnologies
- Cleaning
- Crime scene evaluation
- Dairy processing
- Food Processing
- General maintenance
- Handling small parts
- Janitorial and maintenance work
- Laboratory analysis
- Laboratory testing
- Light duty maintenance and clean up
- Light assembly tasks
- Maintaining equipment & instruments
- Pharmaceutical manufacturing
- Preparing pharmaceutical products
- Protection from blood and other bodily fluids
- Sanitation and dish washing
- Under glove for meat/fish processing

### TECHNICAL DATA SHEET

#### PRODUCT INFORMATION

	Diamond Grip Plus™
Material	Latex
Color	Natural
Glove Design	Ambidextrous, Powder-Free, Chlorinated, Textured Grip
Cuff	Beaded Cuff
Manufacturing/QMS Audit Standards	ISO 13485
Regulatory/Standards Compliance	PPED 89/686/EEC Category III Complex Design, EN 455 1,2 and 3, EN374, EN374-2 level 2 AQL 1.5, ASTM D3578, EN420 (sizing for special purposes), REACH compliant, TGA, US QSR/FDA 510(k), FOOD MIGRATION
Packaging	100 gloves per dispenser, 10 dispensers per shipper carton, 1000 gloves per shipper carton.
Storage	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Country of Origin	Malaysia

#### PHYSICAL PROPERTIES

	TYPICAL VALUES					TESTING METHOD
	5.5 - 6 XS	6.5 - 7 S	7.5 - 8 M	8.5 - 9 L	9.5 - 10 XL	
Length (mm)	245	245	245	245	245	ASTM D3767/EN 420
Average Palm Width (mm)	76	85	95	106	115	
Freedom from Holes (Inspection level I)	1.5 AQL					ASTM D5151/EN 374-2
Palm Thickness Single Wall	mm : 0.13 / mil : 5.1					ASTM D3767/EN 420
	BEFORE AGING		AFTER AGING			
Ultimate Tensile Strength	> 22 Mpa		> 20 Mpa		ASTM D 412-06a	
Elongation at Break (%)	> 700		> 600		ASTM D 412-06a	
Force at break (N)	9		9		EN 455-2	

#### ORDERING INFORMATION

	SIZE	5.5 - 6 XS	6.5 - 7 S	7.5 - 8 M	8.5 - 9 L	9.5 - 10 XL
Diamond Grip Plus™	REORDER NO.	DGP-350-XS	DGP-350-S	DGP-350-M	DGP-350-L	DGP-350-XL

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application.

Warning: This product contains natural rubber latex, which may cause allergic reactions in sensitive individuals.

For additional information visit us at [www.ansell.com](http://www.ansell.com), or call us at

#### Europe, Middle East & Africa Region

Ansell Healthcare Europe NV  
Riverside Business Park  
Blvd International, 55  
1070 Brussels, Belgium  
T: +32 2 528 74 00  
F: +32 2 528 74 01

#### Asia Pacific Region

Ansell Global Trading Center  
(Malaysia) Sdn Bhd  
Prima 6, Prima Avenue  
Block 3512, Jalan Teknokrat 6  
63000 Cyberjaya  
Selangor, Malaysia  
T: +60 3 8310 6688  
F: +60 3 8318 6699

#### North America Region

Ansell Healthcare Products LLC  
111 Wood Avenue South,  
Suite 210  
Iselin, NJ 08830, USA  
T: +1 800 800 0444  
F: +1 800 800 0445

#### Latin America & Caribbean Region

Ansell Commercial México S.A. de C.V.  
Bvld. Bernardo Quintana No. 7001-C,  
Q7001 Torre II.  
Suites 1304, 1305 y 1306.  
Col. Centro Sur, C.P. 76079  
Querétaro, Gro. México  
T: +52 442 248 1544 / 248 3133

#### Australia

Ansell Limited  
Level 3, 678 Victoria Street,  
Richmond, Vic, 3121  
Australia  
T: +61 1800 337 041  
F: +61 1800 803 578

#### Russia

Анселл РУС  
123610 Россия, Москва  
Краснопресненская  
Наб. 12, п.3, оф. 1304А  
Тел. +7 495 258 13 16



EN 374

EN 374



CE PPE 89/686-0493  
MDD 93/42