

SUPRENO[®] SE Durable, Heavy-Duty Coverage Nitrile Exam Gloves

Provides sturdy coverage for demanding applications





Multi-Purpose

FEATURES AND BENEFITS:

Exceptional Protection

The Microflex Supreno® SE offers an advanced barrier protection against contaminants due to its lower acceptable pinhole rates. Its 0.65 Acceptance Quality Limit for pinholes help set a new standard for protection.

Increased Productivity

Non-stick properties are built into the glove formulation to create a resistance to tape or adhesives, and to assist workers with improved efficiency and speed.

Secure and Confident Grip

Textured fingers ensure strong, efficient grip and maximum worker protection.

Effortless Donning and Doffing

Polymer coating assures easy donning to help protect workers in fast-paced environments.

Protection from Type 1 Allergies

Constructed without natural rubber latex materials or powder, to help protect from allergies, skin irritation and dryness.

Industries

- Automotive
- Chemical
- Healthcare
- Laboratory
- Law Enforcement
- Manufacturing

Applications

- Crime scene evaluation
- Patient examinations
- Analytical testing
- Assembly and inspection
- Electronics
- Food processing
- Glass manufacturing
- Laboratory analysis and testing
- Light assembly tasks
- Maintaining equipment & instruments
- Paint and spray shops
- Patient examinations
- Pharmaceutical manufacturing
- Preparing pharmaceutical products
- Printing industry
- Protection from blood and other bodily fluids



SUPRENO[®] SE Durable, Heavy-Duty Coverage Nitrile Exam Gloves

TECHNICAL DATA SHEET:

PRODUCT INFORMATION:

	SUPRENO* SE
Material	Nitrile
Color	Blue
Glove Design	Ambidextrous, Powder-Free, Chlorinated, Polymer Coated, Textured Fingers
Cuff	Beaded Cuff
Audit Standards	ISO 13485
Regulatory Compliance	MDD 93/42/EEC, PPED 89/686/EEC CATEGORY III, COMPLEX DESIGN, EN 455 1,2 AND 3, EN374-3, EN374-2 LEVEL 3 AQL .65, ASTM D6319, EN420 (SIZING FOR SPECIAL PURPOSES), REACH COMPLIANT, TGA, US QSR/FDA 510(K)
Packaging	100 gloves per dispenser, 10 dispensers per shipper carton, 1000 gloves per shipper carton.
Storage	Keep out of direct sunlight; store in a cool and dry place. Keep away from sources of ozone or ignition.
Country of Origin	Malaysia

PHYSICAL PROPERTIES:

		TY	TESTING METHOD				
	5.5 - 6 XS	6.5 - 7 S	7.5 - 8 M	8.5 - 9 ∟	9.5 - 10 XL		
Length (mm)	245	245	245	245	245	ASTM D3767/EN 420	
Average Palm Width (mm)	75	86	98	108	115		
Freedom from Holes (Inspection level I)			ASTM D5151/EN 374-2				
Palm Thickness Single Wall	(mm : 0.11) / (mil : 4.3)					ASTM D3767/EN 420	
	BEF	ORE AGING		AFTER AC	GING		
Ultimate Tensile Strength	;	30 Mpa		30 Mp	a	ASTM D 412-06a	
Elongation at Break (%)	> 500			> 400)	ASTM D 412-06a	

ORDERING INFORMATION:

Force at break (N)

	SIZE	5.5 - 6 XS	6.5 - 7 S	7.5 - 8 M	8.5 - 9 L	9.5 - 10 XL
SUPRENO [®] SE	REORDER NO.	SU-690-XS	SU-690-S	SU-690-M	SU-690-L	SU-690-XL

Neither this document nor any other statement made herein by or on behalf of Ansell should be construed as a warranty of merchantability or that any Ansell product is fit for a particular purpose. Ansell assumes no responsibility for the suitability or adequacy of an end user's selection of gloves for a specific application

For additional information visit us at www.ansell.com, or call us at

Europe, Middle East & Africa

Ansell Healthcare Europe NV Riverside Business Park Blvd International, 55 1070 Brussels, Belgium +32 2 528 74 00 +32 2 528 74 01

Asia Pacific Ansell (Hong Kong) Limited Unit 3102-3103, 31/F Manhattan Place 23 Wang Tai Road Kowloon Bay, Hong Kong +852 2185 0600 +852 2856 2155

North America

Ansell Protective Products Inc.

10

111 Wood Avenue South Suite 210 Iselin, NJ 08830, USA +1 800 800 0444

+1 800 800 0445

Latin America & Caribbean Ansell Brazil Ltda. Rua das Figueiras 474 - 4º Andar Bairro Jardim SP 09080-300 Santo André Brazil CNPJ: 03.496.778/0001-21 +55 11 3356 3100

Australia

11

Ansell Limited Level 3, 678 Victoria Street, Richmond, Vic, 3121 Australia +61 1800 337 041 +61 1800 803 578

Russia Анселл РУС 123610 Россия, Москва Краснопресненская Наб. 12, п.3, оф. 1304A Тел. +7 495 258 13 16





* and ™ are trademarks owned by Ansell Limited or one of its affiliates. © 2014 All Rights Reserved.