

《 COMMODITY DESCRIPTION 》

카탈로그 번호	품목번호	품목 및 규격	단위	수량
Cat. No	Item No	Description	Unit	Quantity
B04-811-259	LO-BG200	Digital Glass Oil Bath	set	1

A. Feature

1. Bath가 투명한 유리로 되어있어 내부 관찰이 용이함
2. Magnetic Stirrer와 같이 사용하여 온도 분포도를 향상 시킬 수 있음
3. 0.01 °C 온도 표시
4. 외부센서 및 RS485 PC 통신용 포트 기본 장착
5. Graphic LCD Display, 간편하고 직관적인 Controller 조작
6. Auto-Tuning 기능 및 온도 Calibration 기능 장착
7. 2가지 타입의 타이머 (구동과 동시에 시작 & 목표 온도 도달 후 시작)
8. Automatical Fault Sensing System 기기 스스로 이상 상황을 감시 (특허)
9. 2중화 과온 안전장치, 1차 Controller 과온 Alarm과 2차 과온 차단 회로

B. Specifications

Model	LO-BG200	
Capacity	4.5 L	
Controller	Control	PID Control, Auto-tuning
	Display	GLCD (Graphic LCD)
	Resolution	0.01 °C
	Timer	Run Start / Wait Start (99 hr 59 min)
Temperature	Range	Ambient +5 to +250 °C
	Accuracy	±0.3 °C (at 100 °C)
	Uniformity	±1.52 °C (at 100 °C)
Dimension	Bath	Φ200×h165 mm
	Packing Dimension	w435 × d295 × h290 mm
Weight	Weight	3.0 kg
	Packing Weight	3.6 kg
Electric Supply	Power	1 Phase, 220 VAC, 50/60 Hz
	Max Consumption	550 W (2.5 A)
	Power Line	Standard Plug
Safety	Over Temp. Limiter & Current Protection, Automatical Fault Sensing System	


C. Accessory

1. B04-91-255 / GB045 / Spare Glass Bath (1ea 기본포함)
2. B04-01-085 / LI-USB / RS485 to USB Communication Cable (옵션)

3. H03-03-205 / SL-HSKIT / Support Kit (옵션)
4. B04-91-471 / ET101 / External Sensor (옵션)

D. Remark

1. 사용자의 지정장소에 설치 후 1년 동안 무상보증.
2. PL Insurance (제조물배상책임보험)
3. Validation Service (IQ, OQ, PQ)

 www.lklab.com	(주)엘케이랩코리아	TEL : 031-573-4952 FAX : 031-571-4958	E-mail : info@lklab.com
---	------------	--	----------------------------