

## 《 COMMODITY DESCRIPTION 》

카탈로그 번호	품목번호	품목 및 규격	단위	수량
Cat. No	Item No	Description	Unit	Quantity
B04-01-230	CTB40	Cold Trap Bath	set	1

### A. Feature

1. 진공라인에서 발생하는 증기 및 유해성 물질 응축용
2. 온도 조절 가능
3. 설치가 용이한 2개의 Glass Trap 사용(기본 포함)
4. 컴팩트 저소음형 동결 트랩
5. 착탈식 전면 커버 구조, 냉동기 유지 관리 용이, 효율 증대
6. Bath 내부 용액을 편리하게 배출시킬 수 있는 Drain 장착
7. Graphic LCD Display, 간편하고 직관적인 Controller 조작
8. 온도 Calibration 기능
9. 2가지 타입의 타이머 (구동과 동시에 시작 & 목표 온도 도달 후 시작)
10. Automatical Fault Sensing System 기기 스스로 이상 상황을 감시 (특허)

### B. Specifications

Model		CTB40
Capacity		4.2 L
Controller	Control	Digital Control
	Resolution	0.1 °C
	Timer	Run Start / Wait Start (99 hr 59 min)
Temperature	Range	Down to -40 °C
	Accuracy	±0.3 °C
	Uniformity	±2.53 °C
Dimension	Interior	Ø210×h220 mm
	External	w340 × d440 × h640 mm
	Packing Dimension	w400 × d470 × h780 mm
Weight	Weight	38 kg
	Packing Weight	41 kg
Cooling	Compressor	1/2 hp
Material	Interior	Stainless Steel 304
	External	Powder Coated Steel
	Cover	PMMA
	Gasket	Silicon Rubber
Electric Supply	Power	1 Phase, 220 VAC, 60 Hz
	Max Consumption	400 W (1.9 A)


	Power Line	Standard Plug
Other	Drain	1/4" Hose Nipple Male
	Caster	Leveling Caster (10 mm)
Safety		Over Current Protection, Automatical Fault Sensing System

**C. Accessory**

1. B04-01-A30 / GCT01 / Glass Trap (2ea 기본 포함)
2. B04-01-A40 / GCT02 / Spare Rubber Gasket (옵션)
3. B04-01-A50 / GCT03 / Viton Gasket (옵션)

**D. Remark**

1. 사용자의 지정장소에 설치 후 1년 동안 무상보증.
2. PL Insurance (제조물배상책임보험)
3. Validation Service (IQ, OQ, PQ)
4. CE Mark
5. KC Mark

 <p>www.lklab.com</p>	<p>(주)엘케이랩코리아</p>	<p>TEL : 031-573-4952 FAX : 031-571-4958</p>	<p>E-mail : info@lklab.com</p>
---	-------------------	--	------------------------------------